

Anti-Lama4 Rabbit Polyclonal Antibody For Wb, Elisa - Q16363

Item Number: CM0014269

Introduction

LAMA4 Rabbit Polyclonal Antibody specifically detects human laminin subunit alpha-4 (203kDa) in WB and ELISA applications. Suitable for studying laminin-mediated cell adhesion, embryonic development, and basement membrane biology.

[Learn More](#)

Parameter	Value
Product Name	LAMA4 Rabbit pAb
Remarks/Alias	LAMA3; CMD1JJ; LAMA4*-1; LAMA4
Species	Human
GeneID (human)	3910
GeneID	3910
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16363
Protein Weight	203kDa
Shipping	Ice bag