

Anti-Lamb3 Rabbit Polyclonal Antibody For Wb, Elisa - Q13751

Item Number: CM0013247

Introduction

LAMB3 rabbit polyclonal antibody (pAb) validated for WB and ELISA with human and rat reactivity. Target laminin subunit beta-3, a basement membrane glycoprotein involved in cell adhesion. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	LAMB3 Rabbit pAb
Remarks/Alias	AI1A; LAM5; JEB1A; JEB1B; LAMNB1; BM600-125KDA; LAMB3
Species	Human
Gene ID	3914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 762-1004 of human LAMB3 (NP_000219.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13751
Protein Weight	130kDa
Shipping	Ice bag