

Anti-Lambda Light Chain Monoclonal Antibody For Wb, Ihc-P, Elisa - P01701

Item Number: CM0003015

Introduction

Rabbit monoclonal antibody targeting human IGLV1-51 (Lambda Light Chain variable region). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for immunoglobulin research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Lambda Light Chain Rabbit mAb
Remarks/Alias	V1-19; IGLV151; Lambda Light chain
Species	Human
Gene ID (Human)	28820
Gene ID	28820
Immunogen	Recombinant protein (aa 1-117 of human Lambda Light chain, SWISS: P01701)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01701
Protein Weight	12kDa
Shipping	Ice bag