


Anti-Lamc1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11047

Item Number: CM0006376



Introduction

LAMC1 rabbit monoclonal antibody validated for WB, IHC-P, ELISA with reactivity in human, mouse, rat. Laminin subunit gamma-1 involved in cell adhesion and basement membrane formation. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	LAMC1 Rabbit mAb
Remarks/Alias	LAMB2
Species	Human
GeneID (human)	3915
GeneID	3915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1030-1275 of human LAMC1 (NP_002284.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11047
Protein Weight	178kDa
Shipping	Ice bag