


Anti-Lamc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, And Elisa - P11047

Item Number: CM0013635



Introduction

High-quality anti-LAMC1 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA applications. Reacts with human, mouse, and rat laminin gamma-1 subunit. Ideal for extracellular matrix and basement membrane research.

[Learn More](#)

Parameter	Value
Product Name	LAMC1 Rabbit pAb
Remarks/Alias	LAMB2; LAMC1
Species	Human
GeneID (Human)	3915
GeneID	3915
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11047
Protein Weight	178kDa
Shipping	Ice bag