

Anti-Lamc2 Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q13753

Item Number: CM0020752

Introduction

Anti-LAMC2 rabbit polyclonal antibody for Western blot, IF/ICC, IF-P, and ELISA; reacts with human, mouse, rat; targets laminin subunit gamma-2, a basement membrane component involved in cell adhesion and scattering.

[Learn More](#)

Parameter	Value
Product Name	LAMC2 Rabbit pAb
Remarks/Alias	B2T; CSF; EBR2; BM600; EBR2A; JEB3A; JEB3B; LAMB2T; LAMNB2; LAMC2
Species	Human
GeneID	3918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 746-860 of human LAMC2 (NP_005553.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13753
Protein Weight	131kDa
Shipping	Ice bag