

Anti-Lamin A/C Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P02545

Item Number: CM0002316



Introduction

[KO Validated] Rabbit monoclonal antibody targeting Lamin A/C. Validated for Western blot, IHC, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, rat. Ideal for studies on nuclear lamina, laminopathies, and cellular senescence.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Lamin A/C Rabbit mAb
Remarks/Alias	FPL; IDC; LFP; CDDC; EMD2; FPLD; HGPS; LDP1; LMN1; LMNC; MADA; PRO1; CDCD1; CMD1A; FPLD2; LMNL1; CMT2B1; LGMD1B; /C
Species	Human
GeneID (Human)	4000
GeneID	4000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 403-572 of human Lamin A/C (NP_733821.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02545
Protein Weight	74kDa
Shipping	Ice bag