


Anti-Lamin B1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - P20700

Item Number: CM0023094



Introduction

Lamin B1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP applications. Detects Human, Mouse, Rat Lamin B1. Ideal for studying nuclear lamina architecture and chromatin organization.

[Learn More](#)

Parameter	Value
Product Name	Lamin B1 Rabbit pAb
Remarks/Alias	LMN; ADLD; LMN2; LMNB; MCPH26; Lamin B1
Species	Human
Gene ID (Human)	4001
Gene ID	4001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 397-586 of human Lamin B1 (NP_005564.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P20700
Protein Weight	66kDa
Shipping	Ice bag