


Anti-Lamin B2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q03252

Item Number: CM0025878



Introduction

Anti-Lamin B2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat Lamin B2 (70 kDa). Ideal for nuclear lamina and cell cycle research.

[Learn More](#)

Parameter Name	Value
Product Name	Lamin B2 Rabbit pAb
Remarks/Alias	EPM9; LMN2; LAMB2; MCPH27; Lamin B2
Species	Human
GeneID (Human)	84823
GeneID	84823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-620 of human Lamin B2 (NP_116126.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03252
Protein Weight	70kDa
Shipping	Ice bag