


Anti-Laminin Beta 1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P07942

Item Number: CM0011110



Introduction

Rabbit monoclonal antibody against human Laminin beta 1, applicable for Western blot, immunocytochemistry/immunofluorescence, and ELISA. Reacts with human, mouse, and rat. Validated for research use in extracellular matrix and basement membrane studies.

[Learn More](#)

Parameter	Value
Product Name	Laminin beta 1 Rabbit mAb
Remarks/Alias	CLM; LIS5; Laminin beta 1
Species	Human
Gene ID (Human)	3912
Gene ID	3912
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1687-1786 of human Laminin beta 1 (P07942).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07942
Protein Weight	198kDa
Shipping	Ice bag