

Anti-Laminin Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q9Y6N6 / P25391 / P24043 / P07942 / P55268 / P11047

Item Number: CM0025620

Introduction

High-quality rabbit polyclonal antibody targeting Laminin subunit gamma-3 and other laminin chains. Validated for WB, IF/ICC, IF-P, IHC-P, ELISA, with reactivity in Human, Mouse, Rat. Ideal for extracellular matrix and basement membrane investigations.

[Learn More](#)

Parameter	Value
Product Name	Laminin Rabbit pAb
Remarks/Alias	ARP; Arp1; ARL18
Species	Human
Gene ID (Human)	10319/284217/3908/3912/3913/3915/3918
Gene ID	10319/284217/3908/3912/3913/3915/3918
Immunogen	Recombinant protein of human Laminin
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6N6/P25391/P24043/P07942/P55268/P11047
Protein Weight	171kDa/337kDa/344kDa/198kDa/196kDa/178kDa
Shipping	Ice bag