


Anti-Lamp1 Polyclonal Antibody For Wb And Elisa - P11279

Item Number: CM0028877



Introduction

Rabbit polyclonal antibody against human LAMP1 (CD107a) for Western blot and ELISA applications. Detects a 45 kDa lysosomal membrane glycoprotein critical for autophagy, NK cell cytotoxicity, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	LAMP1 Rabbit pAb
Remarks/Alias	LAMPA; CD107a; LGP120; LAMP1
Species	Human
Gene ID (Human)	3916
Gene ID	3916
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human LAMP1 (NP_005552.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P11279
Protein Weight	45 kDa
Shipping	Ice bag