


Anti-Lamp1 Polyclonal Antibody (Kd Validated) For Wb, If/Icc, Ihc-P, Elisa - P11279

Item Number: CM0019576



Introduction

Rabbit polyclonal LAMP1 antibody validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Targets lysosomal membrane glycoprotein involved in autophagy, lysosomal pH regulation, and NK cell cytotoxicity. Ideal for cell biology and immuno-oncology research. - P11279

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LAMP1 Rabbit pAb
Remarks/Alias	LAMPA; CD107a; LGP120; LAMP1
Species	Human
GeneID (Human)	3916
GeneID	3916
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human LAMP1 (NP_005552.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11279
Protein Weight	45 kDa
Shipping	Ice bag