


Anti-Lamp1/Cd107A Monoclonal Antibody For Flow Cytometry - P11279

Item Number: CM0008531



Introduction

Rabbit anti-Human LAMP1/CD107a monoclonal antibody validated for intracellular flow cytometry. Targets lysosomal membrane glycoprotein involved in autophagy and NK cell cytotoxicity. Recombinant immunogen (aa 29-382), UniProt P11279, ~38-45 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP1/CD107a mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
GeneID (Human)	3916
GeneID	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag