


Anti-Lamp1/Cd107A Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P11279

Item Number: CM0006711



Introduction

Recombinant rabbit monoclonal antibody targeting human LAMP1/CD107a, validated for WB, IHC-P, IF/ICC, FC, and ELISA. Ideal for lysosome, autophagy, and NK cell cytotoxicity research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LAMP1/CD107a Rabbit mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
GeneID (Human)	3916
GeneID	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag