


Anti-Lamp1/Cd107A Rabbit Polymab® For Wb, Ihc-P, If/Icc, Elisa - P11438/P11279

Item Number: CM0006795



Introduction

KD-validated rabbit polymorphic monoclonal antibody against human/mouse LAMP1/CD107a for WB, IHC-P, IF/ICC, and ELISA; ideal for lysosomal biogenesis and NK cell cytotoxicity research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LAMP1/CD107a Rabbit PolymAb®
Remarks/Alias	P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120
Species	Human/Mouse
GeneID (Human)	3916
GeneID (Mouse)	16783
Immunogen	Recombinant protein of human/mouse CD107a
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11438/P11279
Protein Weight	43kDa
Shipping	Ice bag