

Anti-Lamp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P13473

Item Number: CM0022922



Introduction

Rabbit polyclonal antibody against human LAMP2, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for lysosome biogenesis, autophagy, and MHC class II presentation research.

[Learn More](#)

Parameter	Value
Product Name	LAMP2 Rabbit pAb
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110; LAMP2
Species	Human
GeneID (Human)	3920
GeneID	3920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13473
Protein Weight	45kDa
Shipping	Ice bag