

Anti-Lamp2 Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - P13473

Item Number: CM0015401



Introduction

Rabbit polyclonal antibody targeting LAMP2 (CD107b), validated for Western blot, IHC-P, and ELISA. Cross-reacts with human and mouse. Recombinant immunogen (aa 29-375). Ideal for autophagy, lysosome, and chaperone-mediated autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] LAMP2 Rabbit pAb
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110; P2
Species	Human
Gene ID	3920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13473
Protein Weight	45kDa
Shipping	Ice bag