

Anti-Lamp3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uqv4

Item Number: CM0024706



Introduction

LAMP3 rabbit polyclonal antibody validated for WB, IHC-P, ELISA with reactivity across human, mouse, rat. Suitable for autophagy, dendritic cell, and viral infection research. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	LAMP3 Rabbit pAb
Remarks/Alias	LAMP, CD208, DCLAMP, LAMP-3, TSC403, DC LAMP, DC-LAMP
Species	Human
Gene ID	27074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-380 of human LAMP3 (NP_055213.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQV4
Protein Weight	44kDa
Shipping	Ice bag