

Anti-Lamtor1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6laa8

Item Number: CM0021948



Introduction

Rabbit polyclonal antibody targeting LAMTOR1, validated in WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Suitable for mTOR signaling, lysosome biology, and amino acid sensing research.

[Learn More](#)

Parameter	Value
Product Name	LAMTOR1 Rabbit pAb
Remarks/Alias	p18; PDRO; C11orf59; p27RF-Rho; Ragulator1; LAMTOR1
Species	Human
GeneID (Human)	55004
GeneID	55004
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human LAMTOR1 (NP_060377.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6IAA8
Protein Weight	18kDa
Shipping	Ice bag