


Anti-Lap2 Rabbit Polyclonal Antibody For Wb, Elisa - P42166/P42167

Item Number: CM0022093



Introduction

High-quality rabbit polyclonal antibody against human LAP2 (thymopoietin, TMPO). Validated for Western blot and ELISA. Cross-reacts with mouse and rat. Supports research on nuclear architecture and T-cell development.

[Learn More](#)

Parameter	Value
Product Name	LAP2 Rabbit pAb
Remarks/Alias	TMPO; CMD1T; LAP2; LEMD4; PRO0868; TP; thymopoietin
Species	Human
Gene ID (Human)	7112
Gene ID	7112
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human LAP2 (NP_001027454.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P42166/P42167
Protein Weight	38kDa/50kDa/75kDa
Shipping	Ice bag