

Anti-Lap3 (Kd Validated) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P28838

Item Number: CM0010893

Introduction

KD-validated anti-LAP3 rabbit monoclonal antibody for Western blotting, IHC-P, and ELISA. Detects human, mouse, rat cytosolic aminopeptidase LAP3, involved in glutathione metabolism. Ideal for knockdown validation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LAP3 Rabbit mAb
Remarks/Alias	LAP; PEPS; LAPEP; HEL-S-106
Species	Human
GeneID (human)	51056
GeneID	51056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-180 of human LAP3 (NP_056991.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28838
Protein Weight	56kDa
Shipping	Ice bag