

Anti-Lap3 Rabbit Polyclonal Antibody For Wb, Elisa - P28838

Item Number: CM0025701

Introduction

High-quality rabbit polyclonal antibody against human LAP3 (Cytosol aminopeptidase) for Western blot and ELISA. Reacts with human, mouse, rat. Ideal for studying glutathione metabolism and cell redox status. 56kDa.

[Learn More](#)

Parameter	Value
Product Name	LAP3 Rabbit pAb
Remarks/Alias	LAP; PEPS; LAPEP; HEL-S-106; LAP3
Species	Human
GeneID (human)	51056
GeneID	51056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-519 of human LAP3 (NP_056991.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28838
Protein Weight	56kDa
Shipping	Ice bag