

Anti-Lap3 Rabbit Polyclonal Antibody - P28838

Item Number: CM0020365

Introduction

Rabbit polyclonal anti-human LAP3 antibody. Tested for WB, IF-P, ELISA; cross-reacts with mouse, rat. Investigate cytosolic aminopeptidase function in glutathione metabolism.

[Learn More](#)

Parameter	Value
Product Name	LAP3 Rabbit pAb
Remarks/Alias	LAP3; HEL-S-106; LAP; LAPEP; PEP5; cytosol aminopeptidase
Species	Human
Gene ID	51056
Immunogen	Recombinant protein of human LAP3
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28838
Protein Weight	52kDa/56kDa
Shipping	Ice bag