

Anti-Laptm4B Rabbit Polyclonal Antibody For Wb, Elisa - Q86Vi4

Item Number: CM0015528

Introduction

Rabbit polyclonal antibody targeting LAPTM4B (lysosomal-associated transmembrane protein 4B) for WB and ELISA; detects human, mouse, and rat LAPTM4B, involved in lysosomal function and mTORC1 signaling.

[Learn More](#)

Parameter	Specification
Product Name	LAPTM4B Rabbit pAb
Remarks/Alias	LC27, LAPTM4beta, LAPTM4B
Species	Human
Gene ID (Human)	55353
Gene ID	55353
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 218-317 of human LAPTM4B (NP_060877.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86VI4
Protein Weight	35kDa
Shipping	Ice bag