

Anti-Lar Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P10586

Item Number: CM0012069



Introduction

Human anti-LAR rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Detects ~213 kDa receptor-type tyrosine phosphatase F (PTPRF) involved in cell adhesion and EPHA2 signaling.

[Learn More](#)

Parameter	Value
Product Name	LAR Rabbit mAb
Remarks/Alias	LAR; BNAH2
Species	Human
GeneID (Human)	5792
GeneID	5792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1016-1250 of human PTPRF (NP_002831.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P10586
Protein Weight	213kDa
Shipping	Ice bag