


Anti-Larp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Pkg0

Item Number: CM0026201



Introduction

Rabbit polyclonal antibody against human LARP1, validated for WB and ELISA. LARP1 is an RNA-binding protein that regulates translation of 5'TOP mRNAs downstream of mTORC1, influencing cell growth and proliferation. Ideal for translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	LARP1 Rabbit pAb
Remarks/Alias	LARP; Lar1; Lhp1; LARP1
Species	Human
GeneID (Human)	23367
GeneID	23367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 106-205 of human LARP1 (Q6PKG0).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6PKG0
Protein Weight	124kDa
Shipping	Ice bag