


Anti-Larp4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q71Rc2

Item Number: CM0030825



Introduction

Rabbit polyclonal antibody against human LARP4. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Immunogen: recombinant fusion protein covering amino acids 1-120. Suitable for mRNA translation and cell morphology studies.

[Learn More](#)

Parameter	Value
Product Name	LARP4 Rabbit pAb
Remarks/Alias	PP13296; LARP4
Species	Human
GeneID (Human)	113251
GeneID	113251
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human LARP4 (NP_443111.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q71RC2
Protein Weight	81kDa
Shipping	Ice bag