


Anti-Larp5 Polyclonal Antibody For Wb, Elisa - Q92615

Item Number: CM0014710



Introduction

Rabbit polyclonal antibody against human LARP5 (La-related protein 4B), immunogen aa 519-738. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying mRNA translation regulation.

[Learn More](#)

Parameter	Value
Product Name	LARP5 Rabbit pAb
Remarks/Alias	LARP5; KIAA0217
Species	Human
GeneID (Human)	23185
GeneID	23185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 519-738 of human LARP5 (NP_055970.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92615
Protein Weight	81kDa
Shipping	Ice bag