


Anti-Larp6 Polyclonal Antibody For Wb, Elisa - Q9Brs8

Item Number: CM0021418



Introduction

Rabbit polyclonal antibody against human LARP6, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for collagen translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	LARP6 Rabbit pAb
Remarks/Alias	ACHN; LARP6
Species	Human
Gene ID (Human)	55323
Gene ID	55323
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human LARP6 (NP_060827.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRS8
Protein Weight	55kDa
Shipping	Ice bag