

Anti-Larp7 Rabbit Polyclonal Antibody For Wb, Elisa - Q4G0J3

Item Number: CM0015258



Introduction

Rabbit polyclonal antibody against LARP7 (hLARP7/PIP7S). Validated for WB and ELISA; reacts with human and mouse. Ideal for studying transcription elongation regulation via the 7SK RNP complex. - Q4G0J3

[Learn More](#)

Parameter	Value
Product Name	LARP7 Rabbit pAb
Remarks/Alias	ALAZS; PIP7S; hLARP7; HDCMA18P; LARP7
Species	Human
Gene ID	51574
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-380 of human LARP7 (NP_057732.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q4G0J3
Protein Weight	67kDa
Shipping	Ice bag