

Anti-Lars Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2J5

Item Number: CM0013910



Introduction

Rabbit polyclonal antibody against human LARS (Leucyl-tRNA synthetase, cytoplasmic).

Validated for WB, IF/ICC, ELISA with human and mouse reactivity. Ideal for protein biosynthesis and aminoacyl-tRNA synthetase research.

[Learn More](#)

Parameter	Value
Product Name	LARS Rabbit pAb
Remarks/Alias	LRS; LARS; LEUS; LFIS; ILFS1; LEURS; PIG44; RNTLS; HSPC192; hr025CI
Species	Human
Gene ID (Human)	51520
Gene ID	51520
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P2J5
Protein Weight	134kDa
Shipping	Ice bag