

Anti-Lasp1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q14847

Item Number: CM0029312



Introduction

Rabbit polyclonal antibody targeting human LASP1 (MLN50), validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Detects a 30 kDa protein involved in actin dynamics and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	LASP1 Rabbit pAb
Remarks/Alias	MLN50; Lasp-1; LASP1
Species	Human
GeneID (Human)	3927
GeneID	3927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-205 of human LASP1 (NP_006139.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14847
Protein Weight	30kDa
Shipping	Ice bag