

Anti-Lat Rabbit Polyclonal Antibody For Wb, Elisa - 043561

Item Number: CM0027220

Introduction

High-quality rabbit polyclonal antibody against human LAT (Linker for activation of T-cells family member 1). Validated for Western blot (WB) and ELISA. Ideal for T-cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	LAT Rabbit pAb
Remarks/Alias	LAT1; pp36; IMD52; LAT
Species	Human
GeneID	27040
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 95-262 of human LAT (NP_055202.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43561
Protein Weight	28kDa
Shipping	Ice bag