

Anti-Lat1/Slc7A5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q01650

Item Number: CM0028993



Introduction

Anti-LAT1/SLC7A5 rabbit polyclonal antibody for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets large neutral amino acids transporter small subunit 1; ideal for transporter and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	LAT1/SLC7A5 Rabbit pAb
Remarks/Alias	E16; CD98; LAT1; 4F2LC; MPE16; D16S469E; LAT1/SLC7A5
Species	Human
Gene ID (Human)	8140
Gene ID	8140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human LAT1/SLC7A5 (NP_003477.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01650
Protein Weight	55kDa
Shipping	Ice bag