


Anti-Lats1 Rabbit Monoclonal Antibody For Wb And Elisa - O95835

Item Number: CM0004387



Introduction

High-quality rabbit monoclonal antibody against human LATS1. Validated for WB and ELISA. Ideal for research on Hippo signaling, tumor suppression, and cell cycle regulation. - O95835

[Learn More](#)

Parameter	Value
Product Name	LATS1 Rabbit mAb
Remarks/Alias	wt; WARTS; LATS1
Species	Human
GeneID (Human)	9113
GeneID	9113
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LATS1 (NP_004681.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95835
Protein Weight	127kDa
Shipping	Ice bag