


Anti-Lats1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95835

Item Number: CM0023072



Introduction

Rabbit polyclonal antibody against human LATS1 (O95835). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for Hippo pathway and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	LATS1 Rabbit pAb
Remarks/Alias	wt; WARTS; LATS1
Species	Human
GeneID (Human)	9113
GeneID	9113
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human LATS1 (NP_004681.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95835
Protein Weight	127kDa
Shipping	Ice bag