

Anti-Lats2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrm7

Item Number: CM0020055



Introduction

High-quality rabbit polyclonal antibody against LATS2, a key kinase in the Hippo tumor suppressor pathway. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for oncology and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	LATS2 Rabbit pAb
Remarks/Alias	KPM; LATS2
Species	Human
GeneID (Human)	26524
GeneID	26524
Immunogen	Recombinant Protein; A synthetic peptide corresponding to a sequence within amino acids 70-145 of human LATS2 (NP_055387.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRM7
Protein Weight	120kDa
Shipping	Ice bag