

Anti-Lbp Polyclonal Antibody For Wb, If/Icc, Elisa - P18428

Item Number: CM0022125

Introduction

LBP Rabbit pAb (SWISS P18428) is a polyclonal antibody validated for WB, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Binds bacterial lipopolysaccharides; ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	LBP Rabbit pAb
Remarks/Alias	BPIFD2; LBP
Species	Human
Gene ID (Human)	3929
Gene ID	3929
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human LBP (NP_004130.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18428
Protein Weight	53kDa
Shipping	Ice bag