


Anti-Lbr Polyclonal Antibody For Wb, If/Icc, Elisa - Q14739

Item Number: CM0027051



Introduction

Rabbit polyclonal antibody targeting human LBR (lamin B receptor), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for cholesterol metabolism and nuclear envelope studies.

[Learn More](#)

Parameter	Value
Product Name	LBR Rabbit pAb
Remarks/Alias	PHA; C14SR; LMN2R; PHASK; TDRD18; DHCR14B; LBR
Species	Human
GeneID (Human)	3930
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14739
Protein Weight	71kDa
Shipping	Ice bag