


Anti-Lc3A/Lc3B Rabbit Polyclonal Antibody For Wb - Q9Gzq8/Q9H492

Item Number: CM0030713



Introduction

Rabbit polyclonal antibody targeting LC3A/LC3B, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Useful for studying autophagy and mitophagy. Q9GZQ8/Q9H492.

[Learn More](#)

Specification	Value
Product Name	LC3A/LC3B Rabbit pAb
Remarks/Alias	MAP1LC3B; ATG8F; LC3B; MAP1A/1BLC3; MAP1LC3B-a; MAP1LC3A; ATG8E; LC3; LC3A; MAP1ALC3; MAP1BLC3; LC3A/LC3B
Species	Human
GeneID (Human)	84557/81631
GeneID	84557/81631
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-60 of human LC3A (NP_115903.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8/Q9H492
Protein Weight	14kDa
Shipping	Ice bag