


Anti-Lc3A/Lc3B Rabbit Polymab® For Wb, Ihc-P, If/Icc, Elisa - Q9Gzq8 / Q9H492

Item Number: CM0011757



Introduction

Rabbit monoclonal antibody targeting human LC3A/LC3B for autophagy research. Validated in WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for autophagosome studies.

[Learn More](#)

Parameter	Value
Product Name	LC3A/LC3B Rabbit PolymAb®
Remarks/Alias	LC3; LC3A; ATG8E; MAP1ALC3; MAP1BLC3; LC3A/B; LC3B
Species	Human
Gene ID (Human)	84557/81631
Gene ID	84557/81631
Immunogen	A synthesized peptide derived from human LC3A/LC3B
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8/Q9H492
Protein Weight	14kDa
Shipping	Ice bag