


Anti-Lc3B Monoclonal Antibody For Wb, Elisa - Q9Gzq8

Item Number: CM0002345



Introduction

Mouse monoclonal antibody against human LC3B (MAP1LC3B). Validated for WB and ELISA, reacts with human, mouse, rat. Ideal for autophagy research, detecting the ~15kDa protein.

[Learn More](#)

Parameter	Value
Product Name	LC3B Mouse mAb
Remarks/Alias	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3
Species	Human
GeneID	81631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human LC3B (NP_073729.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8
Protein Weight	15kDa
Shipping	Ice bag