


Anti-Lc3B Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Gzq8

Item Number: CM0003700



Introduction

KD validated rabbit monoclonal antibody targeting LC3B. Suitable for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LC3B Rabbit mAb
Remarks/Alias	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3; 3B
Species	Human
GeneID (Human)	81631
GeneID	81631
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8
Protein Weight	15 kDa
Shipping	Ice bag