


Anti-Lcat Monoclonal Antibody For Wb, Elisa - P04180

Item Number: CM0003438



Introduction

Anti-LCAT rabbit monoclonal antibody targeting human LCAT (P04180). Validated for WB and ELISA with cross-reactivity to mouse and rat. LCAT is essential for cholesterol esterification in lipoproteins. Ideal for cardiovascular and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	LCAT Rabbit mAb
Remarks/Alias	LCAT; lecithin-cholesterol acyltransferase
Species	Human
GeneID (Human)	3931
GeneID	3931
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LCAT (P04180)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04180
Protein Weight	50kDa
Shipping	Ice bag