


## Anti-Lcat Rabbit Polyclonal Antibody For Wb And Elisa - P04180

Item Number: CM0015797



### Introduction

Rabbit polyclonal antibody to human LCAT (lecithin-cholesterol acyltransferase). Validated for WB and ELISA; cross-reacts with mouse. Useful for studying cholesterol metabolism and lipoprotein remodeling.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | LCAT Rabbit pAb                                                                                                  |
| Remarks/Alias    | LCAT                                                                                                             |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 3931                                                                                                             |
| Gene ID          | 3931                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 25-440 of human LCAT (NP_000220.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Mouse                                                                                                            |
| SWISS            | P04180                                                                                                           |
| Protein Weight   | 50kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |