


Anti-Lck Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P06239

Item Number: CM0007681



Introduction

Anti-LCK rabbit monoclonal antibody for WB, IF-P, ELISA. Cross-reacts with human, mouse. Targets tyrosine-protein kinase Lck involved in T-cell signaling. Protein weight 58kDa. UniProt: P06239.

[Learn More](#)

Parameter	Value
Product Name	LCK Rabbit mAb
Remarks/Alias	LSK; YT16; IMD22; p56lck; pp58lck; LCK
Species	Human
GeneID (human)	3932
GeneID	3932
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human LCK (NP_005347.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06239
Protein Weight	58kDa
Shipping	Ice bag