


Anti-Lck Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06239

Item Number: CM0018747



Introduction

High-quality rabbit polyclonal antibody targeting human LCK (p56lck), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying T-cell receptor signaling and lymphocyte activation.

[Learn More](#)

Parameter	Value
Product Name	LCK Rabbit pAb
Remarks/Alias	LSK; YT16; IMD22; p56lck; pp58lck; LCK
Species	Human
Gene ID (Human)	3932
Gene ID	3932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human LCK (NP_005347.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06239
Protein Weight	58kDa
Shipping	Ice bag