


Anti-Lcp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q13094

Item Number: CM0021286



Introduction

Rabbit polyclonal antibody against human LCP2 (SLP-76). Suitable for Western blot and ELISA applications. Cross-reacts with human and mouse. Useful for studying T-cell receptor signaling and immune cell activation.

[Learn More](#)

Parameter	Specification
Product Name	LCP2 Rabbit pAb
Remarks/Alias	IMD81, SLP76, SLP-76, LCP2
Species	Human
Gene ID	3937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 159-533 of human LCP2 (NP_005556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13094
Protein Weight	60kDa
Shipping	Ice bag