


Anti-Ldb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86U70

Item Number: CM0019639



Introduction

Rabbit polyclonal antibody against LDB1 (CLIM2), suited for WB, IHC-P, and ELISA. Recognizes human, mouse, rat LDB1. Ideal for studies of transcriptional regulation and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	LDB1 Rabbit pAb
Remarks/Alias	NLI; CLIM2; LDB-1; CLIM-2; LDB1
Species	Human
GeneID (Human)	8861
GeneID	8861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-50 of human LDB1 (NP_003884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U70
Protein Weight	47kDa
Shipping	Ice bag